

4(4-Aminophenyl)aminopyridine

Art. ID TRC-A637883-100MG
Unit 100 mg
Deliverydetails De Minimis

Description

Application Note: 4(4-Aminophenyl)aminopyridine is a useful intermediate that is used in the preparation of substituted heterocyclic ureas as p38 kinase inhibitors useful in treating cytokine mediated diseases.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4(4-Aminophenyl)aminopyridine	[60172-08-9]				