

3-amino-4-propoxy-Benzoic acid

Art. ID TRC-A699340-250MG
Unit 250 mg
Deliverydetails De Minimis

Description

Application Note: 3-amino-4-propoxy-Benzoic acid is a useful building block used in the synthesis of inhibitors targeting HDAC8 from *Schistosoma mansoni* for treatment of schistosomiasis.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3-amino-4-propoxy-Benzoic acid	[59691-15-5]				