

Copper-alloy CuSn8 - cylinder (40 mm x 30 mm) (BAM-374)

Art. ID ERM-EB374
Unit cylinder
Deliverydetails No Dangerous Good /not restricted

Description

Certified values

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	92,22	%		
	Tin (Sn)	[7440-31-5]	7,6	%		
	Phosphorus (P)	[7723-14-0]	0,17	%		
	Silver (Ag)	[7440-22-4]	0,00121	%		
	Iron (Fe)	[7439-89-6]	0,004	%		
	Manganese (Mn)	[7439-96-5]	0,00043	%		
	Nickel (Ni)	[7440-02-0]	0,00327	%		
	Lead (Pb)	[7439-92-1]	0,00083	%		
	Zinc (Zn)	[7440-66-6]	0,00404	%		