

Spectral Interference Check Solution #3, Varied conc, in 2-5% Nitric Acid

Art. ID ACCU-SICS-03-5
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	30	µg/mL		
	Iron (Fe)	[7439-89-6]	150	µg/mL		
	Nickel (Ni)	[7440-02-0]	20	µg/mL		