

Sulfur in Heavy Mineral Oil (75 cSt), 3000 µg/g in White mineral oil

Art. ID ACCU-SWMO-30X-20ML-PAK
Unit 5 x 20 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Sulfur | [7704-34-9] | 3000 | µg/g | | |