

Sulfur in Light Mineral Oil (20 cSt), 300 µg/g in Light mineral oil

Art. ID ACCU-SWMO-LT-3X-20ML-PAK
Unit 5 x 20 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Sulfur | [7704-34-9] | 10000 | µg/g | | |