

23 components je 100 mg/L in 5% HNO3/0,5% HF

Art. ID LM24-CUS-70277
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

Description

Conservation 6 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	100	mg/L		
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Arsenic (As)	[7440-38-2]	100	mg/L		
	Boron (B)	[7440-42-8]	100	mg/L		
	Barium (Ba)	[7440-39-3]	100	mg/L		
	Beryllium (Be)	[7440-41-7]	100	mg/L		
	Bismuth (Bi)	[7440-69-9]	100	mg/L		
	Cadmium (Cd)	[7440-43-9]	100	mg/L		
	Cobalt (Co)	[7440-48-4]	100	mg/L		
	Chromium (Cr)	[7440-47-3]	100	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Iron (Fe)	[7439-89-6]	100	mg/L		
	Manganese (Mn)	[7439-96-5]	100	mg/L		
	Molybdenum (Mo)	[7439-98-7]	100	mg/L		
	Nickel (Ni)	[7440-02-0]	100	mg/L		
	Lead (Pb)	[7439-92-1]	100	mg/L		
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Selenium (Se)	[7782-49-2]	100	mg/L		
	Tin (Sn)	[7440-31-5]	100	mg/L		
	Titanium (Ti)	[7440-32-6]	100	mg/L		
	Thallium (Tl)	[7440-28-0]	100	mg/L		
	Vanadium (V)	[7440-62-2]	100	mg/L		
	Zinc (Zn)	[7440-66-6]	100	mg/L		