

Custom Solution - 10 components in HNO3 (2%) (CRM, ISO 17034 and ISO 17025)

Art. ID LM24-CUS-70279
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Boron (B)	[7440-42-8]	100	mg/L		
	Calcium (Ca)	[7440-70-2]	10000	mg/L		
	Iron (Fe)	[7439-89-6]	100	mg/L		
	Potassium (K)	[7440-09-7]	5000	mg/L		
	Magnesium (Mg)	[7439-95-4]	2500	mg/L		
	Manganese (Mn)	[7439-96-5]	100	mg/L		
	Sodium (Na)	[7440-23-5]	2500	mg/L		
	Phosphorus (P)	[7723-14-0]	500	mg/L		
	Sulfate (SO42-)		10000	mg/L		