

**Austenitic cast iron - chips (BAS-BCS-CRM 263/3)**

Art. ID                      ECRM-B 451-2  
Unit                         100 g  
Deliverydetails            No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Carbon (C)        | [7440-44-0] | 2,059               | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 2,092               | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 1,079               | %    |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 0,0593              | %    |        |        |
|                  | Sulfur (S)        | [7704-34-9] | 0,0315              | %    |        |        |
|                  | Chromium (Cr)     | [7440-47-3] | 1,097               | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 14,01               | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 6,26                | %    |        |        |