

Cast iron, RM, for solid sample spectrometry, disc 40 x 18 mm

Art. ID CZ-CI-SPL-17-40A
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	3,38	%		
	Manganese (Mn)	[7439-96-5]	0,042	%		
	Silicon (Si)	[7440-21-3]	1,98	%		
	Phosphorus (P)	[7723-14-0]	0,021	%		
	Sulfur (S)	[7704-34-9]	0,0035	%		
	Chromium (Cr)	[7440-47-3]	0,031	%		
	Nickel (Ni)	[7440-02-0]	0,045	%		
	Copper (Cu)	[7440-50-8]	0,01	%		
	Molybdenum (Mo)	[7439-98-7]	0,005	%		
	Magnesium (Mg)	[7439-95-4]	0,007	%		
	Cerium (Ce)	[7440-45-1]	0,012	%		
	Vanadium (V)	[7440-62-2]	0,014	%		
	Titanium (Ti)	[7440-32-6]	0,015	%		
	Aluminium (Al)	[7429-90-5]	0,096	%		
	Tin (Sn)	[7440-31-5]	~0,004	%		
	Boron (B)	[7440-42-8]	0,0008	%		
	Zinc (Zn)	[7440-66-6]	~0,002	%		
	Cobalt (Co)	[7440-48-4]	0,027	%		
	Nitrogen (N)	[7727-37-9]	~0,0063	%		