

Pure copper, cylinder (40 mm x 30 mm)

Art. ID BAM-M385a
Unit cylinder
Deliverydetails No Dangerous Good /not restricted

Description

Certified values

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Matrix	Copper (Cu)	[7440-50-8]	~99,9	%		
	Silver (Ag)	[7440-22-4]	25,4 ± 0,5	mg/kg		
	Aluminium (Al)	[7429-90-5]	13,3 ± 3,2	mg/kg		
	Arsenic (As)	[7440-38-2]	9,4 ± 2	mg/kg		
	Bismuth (Bi)	[7440-69-9]	5,64 ± 0,28	mg/kg		
	Cadmium (Cd)	[7440-43-9]	2,75 ± 0,28	mg/kg		
	Cobalt (Co)	[7440-48-4]	7,4 ± 0,4	mg/kg		
	Chromium (Cr)	[7440-47-3]	10,4 ± 0,8	mg/kg		
	Iron (Fe)	[7439-89-6]	44,2 ± 1,2	mg/kg		
	Manganese (Mn)	[7439-96-5]	9,9 ± 0,9	mg/kg		
	Nickel (Ni)	[7440-02-0]	10,8 ± 0,6	mg/kg		
	Phosphorus (P)	[7723-14-0]	10 ± 1,3	mg/kg		
	Lead (Pb)	[7439-92-1]	10,8 ± 1	mg/kg		
	Antimony (Sb)	[7440-36-0]	14,9 ± 0,8	mg/kg		
	Selenium (Se)	[7782-49-2]	5 ± 0,7	mg/kg		
	Silicon (Si)	[7440-21-3]	7,3 ± 0,8	mg/kg		
	Tin (Sn)	[7440-31-5]	16,1 ± 1,1	mg/kg		
	Tellurium (Te)	[13494-80-9]	8,1 ± 1,2	mg/kg		
	Titanium (Ti)	[7440-32-6]	6,6 ± 1,1	mg/kg		
	Zinc (Zn)	[7440-66-6]	9,2 ± 0,7	mg/kg		