

**PbSb3, disc approx. 40 mm diameter and 40 mm height**

Art. ID BAM-M110  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Description

Certified values

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	0,107 ± 0,008	%		
	Bismuth (Bi)	[7440-69-9]	0,0126 ± 0,0004	%		
	Antimony (Sb)	[7440-36-0]	3,08 ± 0,08	%		
	Selenium (Se)	[7782-49-2]	0,0106 ± 0,0014	%		
	Tin (Sn)	[7440-31-5]	0,131 ± 0,004	%		
	Silver (Ag)	[7440-22-4]	22,6 ± 1,7	mg/kg		
	Copper (Cu)	[7440-50-8]	6,4 ± 0,4	mg/kg		
	Tellurium (Te)	[13494-80-9]	3,8 ± 0,9	mg/kg		