

## **Custom Solution - 13 components 10 µg/mL in Methanol**

Art. ID LM24-CUS-70289  
Unit 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

### Description

CRM, ISO 17034 and ISO 17025, Conservation 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Nitroso-diethanolamin e	[1116-54-7]	10	µg/mL		
	n-Nitrosodimethylamine	[62-75-9]	10	µg/mL		
	N-Nitrosodiethylamine	[55-18-5]	10	µg/mL		
	N-Nitroso-di-n-propylam ine	[621-64-7]	10	µg/mL		
	N-Nitrosodiisopropylami ne	[601-77-4]	10	µg/mL		
	N-Nitroso-di-n-butylami ne	[924-16-3]	10	µg/mL		
	N-Nitroso-di-iso-butyla mine	[997-95-5]	10	µg/mL		
	N-Nitroso-N,N-di(3,5,5- trimethylhexyl)amine	[1207995-62-7]	10	µg/mL		
	4-Nitrosomorpholine	[59-89-2]	10	µg/mL		
	N-Nitrosopiperidine	[100-75-4]	10	µg/mL		
	N-Nitrosodibenzylamine	[5336-53-8]	10	µg/mL		
	N-Nitroso-N-methylanili ne	[614-00-6]	10	µg/mL		
	N-Nitroso-N-ethylanilin e	[612-64-6]	10	µg/mL		