

Trifluoperazine N-beta-D-glucuronide, 100 µg/mL in methanol: water (1:1), analytical standard, ampule of 1 mL (SUPELCO)

Art. ID SAF-T-118-1ML
Unit 1 mL
Deliverydetails Dry ice shipment UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------------------------|---------------|---------------------|------|--------|--------|
| | Trifluoperazine N-beta-D-glucuronide | [165602-90-4] | | | | |