

High Purity Tin, RM, Block 300 g

Art. ID BAS-BCS-RM 192J
Unit block
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Informational Value	Arsenic (As)	[7440-38-2]	<0,0001	%		
Informational Value	Bismuth (Bi)	[7440-69-9]	<0,0001	%		
Informational Value	Carbon (C)	[7440-44-0]	0,001	%		
Informational Value	Cadmium (Cd)	[7440-43-9]	<0,0001	%		
Informational Value	Copper (Cu)	[7440-50-8]	<0,0001	%		
Informational Value	Iron (Fe)	[7439-89-6]	<0,0005	%		
Informational Value	Indium (In)	[7440-74-6]	<0,0001	%		
Informational Value	Nickel (Ni)	[7440-02-0]	<0,000	%		
Informational Value	Lead (Pb)	[7439-92-1]	<0,001	%		
Informational Value	Sulfur (S)	[7704-34-9]	0,0002	%		
Informational Value	Antimony (Sb)	[7440-36-0]	<0,001	%		
Informational Value	Tin (Sn)	[7440-31-5]	99,996	%		
Informational Value	Zinc (Zn)	[7440-66-6]	<0,0001	%		