

Austenitic (Ni-RESIST), Cast Iron, RM, Block 48 mm x 42 mm x 12 mm

Art. ID BAS-NIRM 6/1
Unit block
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Informational Value	Carbon (C)	[7440-44-0]	2,53	%		
Informational Value	Cerium (Ce)	[7440-45-1]	0,006	%		
Informational Value	Chromium (Cr)	[7440-47-3]	1,02	%		
Informational Value	Copper (Cu)	[7440-50-8]	0,11	%		
Informational Value	Manganese (Mn)	[7439-96-5]	4,07	%		
Informational Value	Molybdenum (Mo)	[7439-98-7]	0,51	%		
Informational Value	Nickel (Ni)	[7440-02-0]	26,9	%		
Informational Value	Phosphorus (P)	[7723-14-0]	0,225	%		
Informational Value	Sulfur (S)	[7704-34-9]	0,049	%		
Informational Value	Silicon (Si)	[7440-21-3]	2,68	%		