

3,3',5-Triiodo-L-Thyronine:HCl (T3) (ring-13C6, 99%)100 µg/mL in 0.1N NH3 in Methanol

Art. ID CIL-CLM-7185-C
Unit 1 ML
Deliverydetails UN1230 (6.1) II E2

Description

Chemical purity: 95%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,3',5-Triiodothyronine -(tyrosine phenyl-13C6) hydrochloride	[1217473-60-3]	100	µg/mL		