

**Bisphenol S (13C12, 98%) 100 µg/mL in Methanol**

Art. ID CIL-CLM-9319-1.2  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bisphenol S (13C12,98%)	[80-09-1]/unlab	100	µg/mL		
		elled				