

Fipronil Sulfone (3-Cyano,Pyrazole-3,4,5-13C4,99%, 3-Cyano,5-15N2, 98%) 100 µg/mL in Methanol

Art. ID CIL-CNLM-9643-1.2
Unit 1,2 mL
Deliverydetails UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Fipronil Sulfone (3-Cyano,Pyrazole-3,4,5-13C4,99%, 3-Cyano,5-15N2,98%)	[2446366-94-3]	100	µg/mL		