

**Fipronil Desulfinyl(3-Cyano,Pyrazole-3,4,5- 13C4,99%, 3-Cyano,5-15N2,98%) 100 µg/mL in Methanol**

Art. ID CIL-CNLM-9647-1.2  
Unit 1,2 mL  
Deliverydetails UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Fipronil Desulfinyl (3-Cyano,Pyrazole-3,4,5- 13C4,99%, 3-Cyano,5-15N2,98%)	[2687960-48-9]	100	µg/mL		