

Tefluthrin (D6, 98%) 100 µg/mL in Nonane

Art. ID CIL-DLM-10038-1.2
Unit 1,2 mL
Deliverydetails UN1920 III E1

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------|----------------|---------------------|-------|--------|--------|
| | Tefluthrin (D6,98%) | [2483831-09-8] | 100 | µg/mL | | |