

## **Cypermethrin, Mix Of Stereoisomers (D6, 98%) 100 µg/mL in Nonane**

Art. ID CIL-DLM-10042-1.2  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cypermethrin, Mix Of Stereoisomers (D6,98%)	[82523-65-7]	100	µg/mL		