

Tetradecyl(Tri-N-Butyl)Phosphonium Bromide (D29, 98%) 100 µg/mL in Acetone:Water (75:25)

Art. ID CIL-DLM-9612-1.2
Unit 1,2 mL
Deliverydetails UN1090 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tetradecyl(Tri-N-Butyl) Phosphonium Bromide (D2 9,98%)	[2687960-94-5]	100	µg/mL		