

Mono-Octa PCN Syringe Spike (13C12, 99%) 1000 ng/mL in Nonane

Art. ID CIL-EC-4979-A
Unit 1,2 ML
Deliverydetails UN1920 III E1

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---|---------------|---------------------|-------|--------|--------|
| IUPAC# 9L | 2,5-Dichlorobiphenyl (1 3C12,99%) (IUPAC #9) | [250694-89-4] | 1000 | ng/mL | | |
| IUPAC# 52L | 2,2',5,5'-Tetrachlorobi phenyl (13C12,99%) (IUP AC #52) | [208263-80-3] | 1000 | ng/mL | | |
| IUPAC# 101L | 2,2',4,5,5'-Pentachloro biphenyl (13C12,99%) (I UPAC #101) | [104130-39-4] | 1000 | ng/mL | | |
| IUPAC# 138L | 2,2',3,4,4',5'-Hexachlo robiphenyl (13C12,99%) (IUPAC #138) | [208263-66-5] | 1000 | ng/mL | | |
| IUPAC# 194L | 2,2',3,3',4,4',5,5'-Oct achlorobiphenyl (13C12, 99%) (IUPAC #194) | [208263-74-5] | 1000 | ng/mL | | |