

PCB Mixture in Isooctane PCBs-28/52/101/138/153/180/209 (1:50 Dilution) (unlabelled)

Art. ID CIL-EC-5179/50
Unit 100 µL
Deliverydetails UN1262 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4,4'-Trichlorobiphenyl (PCB 28)	[7012-37-5]	5000	ng/mL		
	2,2',5,5'-Tetrachlorobiphenyl (PCB 52)	[35693-99-3]	5000	ng/mL		
	2,2',4,5,5'-Pentachlorobiphenyl (PCB 101)	[37680-73-2]	5000	ng/mL		
	2,2',3,4,4',5'-Hexachlorobiphenyl (PCB 138)	[35065-28-2]	5000	ng/mL		
	2,2',4,4',5,5'-Hexachlorobiphenyl (PCB 153)	[35065-27-1]	5000	ng/mL		
	2,2',3,4,4',5,5'-Heptachlorobiphenyl (PCB 180)	[35065-29-3]	5000	ng/mL		
	Decachlorobiphenyl (PCB 209)	[2051-24-3]	5000	ng/mL		