

Comprehensive PCB Calibration Solution CS3 ((unlabelled)/13C12, 99%) in Nonane/Isooctane

Art. ID CIL-EC-5518-CS3
Unit 0,2 mL
Deliverydetails UN1993 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
CS3, IUPAC# 1	2-Monochlorobiphenyl (P CB 1)	[2051-60-7]	4	ng/mL		
CS3, IUPAC# 3	4-Monochlorobiphenyl (P CB 3)	[2051-62-9]	4	ng/mL		
CS3, IUPAC# 4	2,2'-Dichlorobiphenyl (P PCB 4)	[13029-08-8]	4	ng/mL		
CS3, IUPAC# 8	2,4'-Dichlorobiphenyl (P PCB 8)	[34883-43-7]	4	ng/mL		
CS3, IUPAC# 9	2,5-Dichlorobiphenyl (P CB 9)	[34883-39-1]	4	ng/mL		
CS3, IUPAC# 10	2,6-Dichlorobiphenyl (P CB 10)	[33146-45-1]	4	ng/mL		
CS3, IUPAC# 11	3,3'-Dichlorobiphenyl (P PCB 11)	[2050-67-1]	4	ng/mL		
CS3, IUPAC# 12	3,4-Dichlorobiphenyl (P CB 12)	[2974-92-7]	4	ng/mL		
CS3, IUPAC# 15	4,4'-Dichlorobiphenyl (P PCB 15)	[2050-68-2]	4	ng/mL		
CS3, IUPAC# 18	2,2',5-Trichlorobipheny I (PCB 18)	[37680-65-2]	2	ng/mL		
CS3, IUPAC# 19	2,2',6-Trichlorobipheny I (PCB 19)	[38444-73-4]	2	ng/mL		
CS3, IUPAC# 28	2,4,4'-Trichlorobipheny I (PCB 28)	[7012-37-5]	2	ng/mL		
CS3, IUPAC# 31	2,4',5-Trichlorobipheny I (PCB 31)	[16606-02-3]	2	ng/mL		
CS3, IUPAC# 33	2',3,4-Trichlorobipheny I (PCB 33)	[38444-86-9]	2	ng/mL		
CS3, IUPAC# 35	3,3',4-Trichlorobipheny	[37680-69-6]	2	ng/mL		

	I (PCB 35)			
CS3, IUPAC# 37	3,4,4'-Trichlorobiphenyl	[38444-90-5]	2	ng/mL
	I (PCB 37)			
CS3, IUPAC# 38	3,4,5-Trichlorobiphenyl	[53555-66-1]	2	ng/mL
	(PCB 38)			
CS3, IUPAC# 44	2,2',3,5'-Tetrachlorobiphenyl (PCB 44)	[41464-39-5]	2	ng/mL
CS3, IUPAC# 49	2,2',4,5'-Tetrachlorobiphenyl (PCB 49)	[41464-40-8]	2	ng/mL
CS3, IUPAC# 52	2,2',5,5'-Tetrachlorobiphenyl (PCB 52)	[35693-99-3]	2	ng/mL
CS3, IUPAC# 54	2,2',6,6'-Tetrachlorobiphenyl (PCB 54)	[15968-05-5]	2	ng/mL
CS3, IUPAC# 57	2,3,3',5'-Tetrachlorobiphenyl (PCB 57)	[70424-67-8]	2	ng/mL
CS3, IUPAC# 66	2,3',4,4'-Tetrachlorobiphenyl (PCB 66)	[32598-10-0]	2	ng/mL
CS3, IUPAC# 70	2,3',4',5'-Tetrachlorobiphenyl (PCB 70)	[32598-11-1]	2	ng/mL
CS3, IUPAC# 74	2,4,4',5'-Tetrachlorobiphenyl (PCB 74)	[32690-93-0]	2	ng/mL
CS3, IUPAC# 77	3,3',4,4'-Tetrachlorobiphenyl (PCB 77)	[32598-13-3]	2	ng/mL
CS3, IUPAC# 78	3,3',4,5'-Tetrachlorobiphenyl (PCB 78)	[70362-49-1]	2	ng/mL
CS3, IUPAC# 79	3,3',4,5'-Tetrachlorobiphenyl (PCB 79)	[41464-48-6]	2	ng/mL
CS3, IUPAC# 81	3,4,4',5'-Tetrachlorobiphenyl (PCB 81)	[70362-50-4]	2	ng/mL
CS3, IUPAC# 87	2,2',3,4,5'-Pentachlorobiphenyl (PCB 87)	[38380-02-8]	2	ng/mL
CS3, IUPAC# 95	2,2',3,5',6'-Pentachlorobiphenyl (PCB 95)	[38379-99-6]	2	ng/mL
CS3, IUPAC# 99	2,2',4,4',5'-Pentachlorobiphenyl (PCB 99)	[38380-01-7]	2	ng/mL
CS3, IUPAC# 101	2,2',4,5,5'-Pentachlorobiphenyl (PCB 101)	[37680-73-2]	2	ng/mL
CS3, IUPAC# 104	2,2',4,6,6'-Pentachlorobiphenyl (PCB 104)	[56558-16-8]	2	ng/mL

	biphenyl (PCB 104)			
CS3, IUPAC# 105	2,3,3',4,4'-Pentachloro biphenyl (PCB 105)	[32598-14-4]	2	ng/mL
CS3, IUPAC# 110	2,3,3',4',6-Pentachloro biphenyl (PCB 110)	[38380-03-9]	2	ng/mL
CS3, IUPAC# 111	2,3,3',5,5'-Pentachloro biphenyl (PCB 111)	[39635-32-0]	2	ng/mL
CS3, IUPAC# 114	2,3,4,4',5-Pentachlorob iphenyl (PCB 114)	[74472-37-0]	2	ng/mL
CS3, IUPAC# 118	2,3',4,4',5-Pentachloro biphenyl (PCB 118)	[31508-00-6]	2	ng/mL
CS3, IUPAC# 123	2',3,4,4',5-Pentachloro biphenyl (PCB 123)	[65510-44-3]	2	ng/mL
CS3, IUPAC# 126	3,3',4,4',5-Pentachloro biphenyl (PCB 126)	[57465-28-8]	2	ng/mL
CS3, IUPAC# 138	2,2',3,4,4',5'-Hexachlo robiphenyl (PCB 138)	[35065-28-2]	2	ng/mL
CS3, IUPAC# 149	2,2',3,4',5',6-Hexachlo robiphenyl (PCB 149)	[38380-04-0]	2	ng/mL
CS3, IUPAC# 153	2,2',4,4',5,5'-Hexachlo robiphenyl (PCB 153)	[35065-27-1]	2	ng/mL
CS3, IUPAC# 155	2,2',4,4',6,6'-Hexachlo robiphenyl (PCB 155)	[33979-03-2]	2	ng/mL
CS3, IUPAC# 156	2,3,3',4,4',5-Hexachlor obiphenyl (PCB 156)	[38380-08-4]	2	ng/mL
CS3, IUPAC# 157	2,3,3',4,4',5'-Hexachlo robiphenyl (PCB 157)	[69782-90-7]	2	ng/mL
CS3, IUPAC# 162	2,3,3',4',5,5'-Hexachlo robiphenyl (PCB 162)	[39635-34-2]	2	ng/mL
CS3, IUPAC# 167	2,3',4,4',5,5'-Hexachlo robiphenyl (PCB 167)	[52663-72-6]	2	ng/mL
CS3, IUPAC# 169	3,3',4,4',5,5'-Hexachlo robiphenyl (PCB 169)	[32774-16-6]	2	ng/mL
CS3, IUPAC# 170	2,2',3,3',4,4',5-Heptac hlorobiphenyl (PCB 170)	[35065-30-6]	2	ng/mL
CS3, IUPAC# 174	2,2',3,3',4,5,6'-Heptac hlorobiphenyl (PCB 174)	[38411-25-5]	2	ng/mL
CS3, IUPAC# 178	2,2',3,3',5,5',6-Heptac	[52663-67-9]	2	ng/mL

	hlorobiphenyl (PCB 178)			
CS3, IUPAC# 180	2,2',3,4,4',5,5'-Heptac	[35065-29-3]	2	ng/mL
	hlorobiphenyl (PCB 180)			
CS3, IUPAC# 187	2,2',3,4',5,5',6-Heptac	[52663-68-0]	2	ng/mL
	hlorobiphenyl (PCB 187)			
CS3, IUPAC# 188	2,2',3,4',5,6,6'-Heptac	[74487-85-7]	2	ng/mL
	hlorobiphenyl (PCB 188)			
CS3, IUPAC# 189	2,3,3',4,4',5,5'-Heptac	[39635-31-9]	2	ng/mL
	hlorobiphenyl (PCB 189)			
CS3, IUPAC# 194	2,2',3,3',4,4',5,5'-Oct	[35694-08-7]	2	ng/mL
	achlorobiphenyl (PCB 194)			
CS3, IUPAC# 195	2,2',3,3',4,4',5,6-Octa	[52663-78-2]	2	ng/mL
	chlorobiphenyl (PCB 195)			
CS3, IUPAC# 200	2,2',3,3',4,5,6,6'-Octa	[52663-73-7]	2	ng/mL
	chlorobiphenyl (IUPAC P CB 200, PCB 199)			
CS3, IUPAC# 202	2,2',3,3',5,5',6,6'-Oct	[2136-99-4]	2	ng/mL
	achlorobiphenyl (PCB 202)			
CS3, IUPAC# 203	2,2',3,4,4',5,5',6-Octa	[52663-76-0]	2	ng/mL
	chlorobiphenyl (PCB 203)			
CS3, IUPAC# 205	2,3,3',4,4',5,5',6-Octa	[74472-53-0]	2	ng/mL
	chlorobiphenyl (PCB 205)			
CS3, IUPAC# 206	2,2',3,3',4,4',5,5',6-N	[40186-72-9]	2	ng/mL
	onachlorobiphenyl (PCB 206)			
CS3, IUPAC# 208	2,2',3,3',4,5,5',6,6'-N	[52663-77-1]	2	ng/mL
	onachlorobiphenyl (PCB 208)			