

PCB Exhibit Native Par Solution (unlabelled)

Art. ID CIL-EC-5534
Unit 1,2 mL

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
IUPAC# 52	2,2',5,5'-Tetrachlorobi phenyl (PCB 52)	[35693-99-3]	200	ng/mL		
IUPAC# 101	2,2',4,5,5'-Pentachloro biphenyl (PCB 101)	[37680-73-2]	200	ng/mL		
IUPAC# 110	2,3,3',4',6-Pentachloro biphenyl (PCB 110)	[38380-03-9]	200	ng/mL		
IUPAC# 153	2,2',4,4',5,5'-Hexachlo robiphenyl (PCB 153)	[35065-27-1]	200	ng/mL		
IUPAC# 187	2,2',3,4',5,5',6-Heptac hlorobiphenyl (PCB 187)	[52663-68-0]	200	ng/mL		
IUPAC# 170	2,2',3,3',4,4',5-Heptac hlorobiphenyl (PCB 170)	[35065-30-6]	200	ng/mL		
IUPAC# 74	2,4,4',5-Tetrachlorobip henyl (PCB 74)	[32690-93-0]	200	ng/mL		
IUPAC# 99	2,2',4,4',5-Pentachloro biphenyl (PCB 99)	[38380-01-7]	200	ng/mL		
IUPAC# 118	2,3',4,4',5-Pentachloro biphenyl (PCB 118)	[31508-00-6]	200	ng/mL		
IUPAC# 105	2,3,3',4,4'-Pentachloro biphenyl (PCB 105)	[32598-14-4]	200	ng/mL		
IUPAC# 156	2,3,3',4,4',5-Hexachlor obiphenyl (PCB 156)	[38380-08-4]	200	ng/mL		
IUPAC# 209	2,2',3,3',4,4',5,5',6,6 '-Decachlorobiphenyl (P CB 209)	[2051-24-3]	200	ng/mL		
IUPAC# 66	2,3',4,4'-Tetrachlorobi phenyl (PCB 66)	[32598-10-0]	200	ng/mL		
IUPAC# 87	2,2',3,4,5'-Pentachloro biphenyl (PCB 87)	[38380-02-8]	200	ng/mL		
IUPAC# 146	2,2',3,4',5,5'-Hexachlo	[51908-16-8]	200	ng/mL		

	robiphenyl (PCB 146)			
IUPAC# 138	2,2',3,4,4',5'-Hexachlo	[35065-28-2]	200	ng/mL
	robiphenyl (PCB 138)			
IUPAC# 180	2,2',3,4,4',5,5'-Heptac	[35065-29-3]	200	ng/mL
	hlorobiphenyl (PCB 180)			