

DL-PCB Native Plus Mixture (unlabelled) in Isooctane

Art. ID CIL-EC-5559
Unit 1,2 mL
Deliverydetails UN1262 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
IUPAC# 77	3,3',4,4'-Tetrachlorobi phenyl (PCB 77)	[32598-13-3]	2000	ng/mL		
IUPAC# 81	3,4,4',5-Tetrachlorobip henyl (PCB 81)	[70362-50-4]	1000	ng/mL		
IUPAC# 105	2,3,3',4,4'-Pentachloro biphenyl (PCB 105)	[32598-14-4]	2000	ng/mL		
IUPAC# 114	2,3,4,4',5-Pentachlorob iphenyl (PCB 114)	[74472-37-0]	1000	ng/mL		
IUPAC# 118	2,3',4,4',5-Pentachloro biphenyl (PCB 118)	[31508-00-6]	5000	ng/mL		
IUPAC# 123	2',3,4,4',5-Pentachloro biphenyl (PCB 123)	[65510-44-3]	1000	ng/mL		
IUPAC# 126	3,3',4,4',5-Pentachloro biphenyl (PCB 126)	[57465-28-8]	1000	ng/mL		
IUPAC# 156	2,3,3',4,4',5-Hexachlor obiphenyl (PCB 156)	[38380-08-4]	2000	ng/mL		
IUPAC# 157	2,3,3',4,4',5'-Hexachlo robiphenyl (PCB 157)	[69782-90-7]	1000	ng/mL		
IUPAC# 167	2,3',4,4',5,5'-Hexachlo robiphenyl (PCB 167)	[52663-72-6]	1000	ng/mL		
IUPAC# 169	3,3',4,4',5,5'-Hexachlo robiphenyl (PCB 169)	[32774-16-6]	1000	ng/mL		
IUPAC# 189	2,3,3',4,4',5,5'-Heptac hlorobiphenyl (PCB 189)	[39635-31-9]	1000	ng/mL		