

Phosphorus Flame Retardant Standard Mixture (D, 98%) 10 µg/mL in Acetonitrile

Art. ID CIL-ES-5529
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tris(2-chloroethyl)phosphate (D12, 98%)	[115-96-8]/unlabelled	10	µg/mL		
	Tris(2-chloroisopropyl)phosphate (D18, 98%)	[13674-84-5]/unlabelled	10	µg/mL		
	Tris(1,3-dichloro-2-propyl)phosphate (D15, 98%)	[13674-87-8]/unlabelled	10	µg/mL		
	Triethyl phosphate (D15, 98%)	[78-40-0]/unlabelled	10	µg/mL		
	Tripropyl phosphate (D21, 98%)	[513-08-6]/unlabelled	10	µg/mL		
	Tributyl phosphate (D27, 98%)	[126-73-8]/unlabelled	10	µg/mL		
	Triphenyl phosphate(D15, 98%)	[1173020-30-8]	10	µg/mL		