

Phosphorus Flame Retardant Standard Mixture (unlabelled) 10 µg/mL in Acetonitrile

Art. ID CIL-ES-5530
Unit 1,2 mL
Deliverydetails UN1648 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tris(2-chloroethy) phosphate		10	µg/mL		
	Tris(2-chloroisopropyl) phosphate		10	µg/mL		
	Tris(1,3-dichloro-2-propyl) phosphate		10	µg/mL		
	Triethyl phosphate	[78-40-0]	10	µg/mL		
	Tripropyl phosphate	[513-08-6]	10	µg/mL		
	Tributyl phosphate	[126-73-8]	10	µg/mL		
	Triphenyl phosphate	[115-86-6]	10	µg/mL		