

## **JECS Labeled Mixture Solution in Water (Ultrapure)**

Art. ID CIL-ES-5535  
Unit 1,2 mL

### Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter                           | CAS number    | Concentration/Value | Unit  | Method | Source |
|------------------|---|---------------|---------------------|-------|--------|--------|
|                  | 8-Hydroxy-2'-deoxyguano<br>sine (15N5, 98%) | [569649-11-2] | 7500                | ng/mL |        |        |
|                  | DI-Cotinine (2',3',4'-1<br>3C3,99%)         | [n/a]         | 150                 | ng/mL |        |        |