

US EPA Method 8276 Toxaphene Composite Stock Standard (unlabelled) in Nonane

Art. ID CIL-ES-5543
Unit 1,2 mL
Deliverydetails UN1920 III E1

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Hx-Sed	[57981-29-0]	1000	ng/mL		
	Hp-Sed	[208049-58-5]	1000	ng/mL		
	Parlar 26 (2-endo,3-exo ,5-endo,6-exo,8,8,10,10 -Octachlorobornane)	[142534-71-2]	1000	ng/mL		
	Parlar 40	[166021-27-8]	1000	ng/mL		
	Parlar 41	[165820-16-6]	1000	ng/mL		
	Parlar 44	[165820-17-7]	1000	ng/mL		
	Parlar 50 (2-endo,3-exo ,5-endo,6-exo,8,8,9,10, 10-Nonachlorobornane)	[66860-80-8]	1000	ng/mL		
	Parlar 62 (2-endo,2-exo ,5-endo,5-exo,8,9,9,10, 10-Nonachlorobornane)	[154159-06-5]	1000	ng/mL		