

Alternate Metals Standard I, Varied conc, in 2-5% Nitric Acid

Art. ID ACCU-WPAM-01-5
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	20	µg/mL		
	Antimony (Sb)	[7440-36-0]	5	µg/mL		
	Beryllium (Be)	[7440-41-7]	5	µg/mL		
	Cobalt (Co)	[7440-48-4]	10	µg/mL		
	Copper (Cu)	[7440-50-8]	10	µg/mL		
	Iron (Fe)	[7439-89-6]	20	µg/mL		
	Manganese (Mn)	[7439-96-5]	10	µg/mL		
	Nickel (Ni)	[7440-02-0]	10	µg/mL		
	Thallium (Tl)	[7440-28-0]	5	µg/mL		
	Vanadium (V)	[7440-62-2]	20	µg/mL		
	Zinc (Zn)	[7440-66-6]	10	µg/mL		