

PFOS/PFOA Calibration Standard CS0.25 ((unlabelled)/13C, 99%) in Methanol

Art. ID CIL-ES-5570-CS0.25
Unit 0,25 mL
Deliverydetails UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
CS0.25, Native	Perfluorooctanoic acid	[335-67-1]	0,5	ng/mL		
CS0.25, Native	Perfluorooctanesulfonat e, sodium salt	[4021-47-0]	0,5	ng/mL		
CS0.25, Extraction	Perfluorooctanoic acid (13C8,99%)		50	ng/mL		
CS0.25, Extraction	Perfluorooctanesulfonat e, sodium salt (13C8,99 %)		50	ng/mL		
CS0.25, Injection	Perfluorononanoic acid (13C9,99%)		50	ng/mL		