

PFOS/PFOA Calibration Standard CS1 (unlabelled)/13C, 99% in Methanol

Art. ID CIL-ES-5570-CS1
Unit 0,25 mL
Deliverydetails UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
CS1, Native	Perfluorooctanoic acid	[335-67-1]	2	ng/mL		
CS1, Native	Perfluorooctanesulfonat e, sodium salt	[4021-47-0]	2	ng/mL		
CS1, Extraction	Perfluorooctanoic acid (13C8,99%)		50	ng/mL		
CS1, Extraction	Perfluorooctanesulfonat e, sodium salt (13C8,99 %)		50	ng/mL		
CS1, Injection	Perfluorononanoic acid (13C9,99%)		50	ng/mL		