

PAH Recovery Standard (D, 98%) 10 ng/mL in Isooctane:Toluene-D8

Art. ID CIL-ES-9464
Unit 10 mL
Deliverydetails UN1262 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acenaphthene (D10,98%)	[15067-26-2]	10	ng/mL		
	p-Terphenyl (D14,98%)	[1718-51-0]	10	ng/mL		
	Perylene (D12,98%)	[1520-96-3]	10	ng/mL		