

2,2',3,4',5,6-Hexachlorobiphenyl (PCB 147) 100 µg/mL in Isooctane ((unlabelled) (Certified Standard))

Art. ID CIL-PCB-147-CS
Unit 1,2 mL
Deliverydetails UN1262 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',3,4',5,6-Hexachlorobiphenyl (PCB 147)	[68194-13-8]	100	µg/mL		