

**2,3,4,4',5,6-Hexachlorobiphenyl (PCB 166) 100 µg/mL in Isooctane ((unlabelled) (Certified Standard)**

Art. ID CIL-PCB-166-CS  
Unit 1,2 mL  
Deliverydetails UN1262 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3,4,4',5,6-Hexachloro biphenyl (PCB 166)	[41411-63-6]	100	µg/mL		