

## **Standard Bullet Replica**

Art. ID	NIST-2460a
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

Standard Reference Material NIST-2460a is a bullet signature standard comprising bullet profile signatures of six Land Engraved Areas (LEAs) from fired bullets. It is replicated from the original SRM 2460 Standard Bullet. NIST-2460a is intended primarily for use as a check standard for forensic laboratories to help verify that the computerized optical equipment for bullet imaging and profiling is operating properly to measure bullet topography signatures. A unit of NIST-2460a consists of a nickel and gold-coated standard bullet replica that is mounted on a plastic holder. Certified for Normalized Cross-Correlation Function Maximum CCFmax and Signature Difference Ds.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Please ask for current certificate						