

Copper base, leaded brass (chill cast), 350, Disc 40x15 mm

Art. ID MBH-31X7835.6D
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0.527	%		
	Arsenic (As)	[7440-38-2]	0.0069	%		
	Bismuth (Bi)	[7440-69-9]	0.0047	%		
	Boron (B)	[7440-42-8]	0.0032	%		
	Cadmium (Cd)	[7440-43-9]	0.0017	%		
	Cobalt (Co)	[7440-48-4]	0.0059	%		
	Copper (Cu)	[7440-50-8]	59.9	%		
	Iron (Fe)	[7439-89-6]	0.118	%		
	Lead (Pb)	[7439-92-1]	1.31	%		
	Nickel (Ni)	[7440-02-0]	0.061	%		
	Phosphorus (P)	[7723-14-0]	0.0380	%		
	Selenium (Se)	[7782-49-2]	~0,001	%		
	Silicon (Si)	[7440-21-3]	0.004	%		
	Silver (Ag)	[7440-22-4]	0.0048	%		
	Tin (Sn)	[7440-31-5]	0.702	%		
	Zinc (Zn)	[7440-66-6]	37.2	%		
	Tellurium (Te)	[13494-80-9]	0.0007	%		