

Copper Base, Typ: Brass - Leaded, CRM, Provisional Analysis, chill cast, Disc ~40 mm

Art. ID MBH-31X7835.3K
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Provisional Composition	Tin (Sn)	[7440-31-5]	0,36	%		
Provisional Composition	Lead (Pb)	[7439-92-1]	1,7	%		
Provisional Composition	Zinc (Zn)	[7440-66-6]	36,6	%		
Provisional Composition	Iron (Fe)	[7439-89-6]	0,49	%		
Provisional Composition	Nickel (Ni)	[7440-02-0]	0,15	%		
Provisional Composition	Aluminium (Al)	[7429-90-5]	0,48	%		
Provisional Composition	Silicon (Si)	[7440-21-3]	0,075	%		
Provisional Composition	Arsenic (As)	[7440-38-2]	0,058	%		
Provisional Composition	Antimony (Sb)	[7440-36-0]	0,061	%		
Provisional Composition	Phosphorus (P)	[7723-14-0]	0,035	%		
Provisional Composition	Bismuth (Bi)	[7440-69-9]	0,029	%		
Provisional Composition	Cobalt (Co)	[7440-48-4]	0,006	%		
Provisional Composition	Silver (Ag)	[7440-22-4]	0,021	%		
Provisional Composition	Tellurium (Te)	[13494-80-9]	0,001	%		
Provisional Composition	Cadmium (Cd)	[7440-43-9]	0,006	%		
Provisional Composition	Copper (Cu)	[7440-50-8]	60	%		