

Sulfur in Mineral Oil Calibration Standard - 11 Standards per set 0-5.0 Wt%, 100 ml each

Art. ID	ASI-SMO11
Unit	11x100 mL
Deliverydetails	No Dangerous Good /not restricted

Description

Sulfur in Mineral Oil Calibration Standards, 11 standards per set, 100 ml each, Concentrations: Blank, 0.0020, 0.0050, 0.0100, 0.0300, 0.0600, 0.1000, 0.5000, 1.0000, 2.5000, 5.0000 Wt %. Sulfur in Petroleum Oil and Petroleum fuels are manufactured and certified gravimetrically under the guidance of ISO 17034

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]	0	%		
	Sulfur (S)	[7704-34-9]	0,002	%		
	Sulfur (S)	[7704-34-9]	0,005	%		
	Sulfur (S)	[7704-34-9]	0,01	%		
	Sulfur (S)	[7704-34-9]	0,03	%		
	Sulfur (S)	[7704-34-9]	0,06	%		
	Sulfur (S)	[7704-34-9]	0,1	%		
	Sulfur (S)	[7704-34-9]	0,5	%		
	Sulfur (S)	[7704-34-9]	1	%		
	Sulfur (S)	[7704-34-9]	2,5	%		
	Sulfur (S)	[7704-34-9]	5	%		