

## **Hydrogen in Carbon CRM**

Art. ID                      AR-6123  
Unit                         25 g  
Deliverydetails            No Dangerous Good /not restricted

### Description

The intended use of this standard is for the verification and quality check of hydrogen in carbon matrix using ASTM methods utilizing inert gas fusion with thermal conductivity or infrared detection instrumentation. Lot: 240216

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Hydrogen (H)	[1333-74-0]	0,419	%		