

Internal Standard Mix, 4,0 mg/mL deuterated PAHs in Methylene chloride (Dichloromethane)

Art. ID ACCU-Z-014J-PAK
Unit 5 x 1 mL
Deliverydetails De Minimis UN1593 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acenaphthene D10	[15067-26-2]	4	mg/mL		
	Chrysene D12	[1719-03-5]	4	mg/mL		
	1,4-Dichlorobenzene D4	[3855-82-1]	4	mg/mL		
	Naphthalene D8	[1146-65-2]	4	mg/mL		
	Perylene D12	[1520-96-3]	4	mg/mL		
	Phenanthrene D10	[1517-22-2]	4	mg/mL		