

Mercury, nitric acid Solution of (acc. to old Ph Eur.)

Art. ID	LM24-OP-2000678
Unit	100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

Description

European Pharmacopoeia / Reagents / These solutions are produced according to the specifications given in chapter 4.1 using Ph Eur reagent quality components and water R meeting the Ph Eur monograph requirements for Purified water. Reagent solutions used in the limit tests for barium, calcium and sulfates are prepared using distilled water R