

Sodium metabisulphite 45.5 g/L Solution

| | |
|-----------------|---|
| Art. ID | LM24-OP-2000700 |
| Unit | 100 mL |
| Deliverydetails | ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) |

Description

European Pharmacopoeia / Reagents / These solutions are produced according to the specifications given in chapter 4.1 using Ph Eur reagent quality components and water R meeting the Ph Eur monograph requirements for Purified water. Reagent solutions used in the limit tests for barium, calcium and sulfates are prepared using distilled water R