

## **TLC performance test Solution**

Art. ID	LM24-OP-2000839
Unit	10 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

### **Description**

European Pharmacopoeia / Reagents / These solutions are produced according to the specifications given in chapter 4.1 using Ph Eur reagent quality components and water R meeting the Ph Eur monograph requirements for Purified water. Reagent solutions used in the limit tests for barium, calcium and sulfates are prepared using distilled water R