

### **Aluminium alloy AlMg3, chips**

Art. ID BAM-300  
Unit 100 g  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	0,14	%		
	Manganese (Mn)	[7439-96-5]	0,018	%		
	Iron (Fe)	[7439-89-6]	0,198	%		
	Copper (Cu)	[7440-50-8]	0,04	%		
	Titanium (Ti)	[7440-32-6]	0,012	%		
	Zinc (Zn)	[7440-66-6]	0,128	%		
	Magnesium (Mg)	[7439-95-4]	2,68	%		
	Lead (Pb)	[7439-92-1]	0,014	%		
	Chromium (Cr)	[7440-47-3]	0,216	%		
	Tin (Sn)	[7440-31-5]	<0,0005	%		
	Aluminium (Al)	[7429-90-5]	rest			