

**Potassium Chloride General and Ion Activity Standard (Dried at 110 °C)**

Art. ID	NIST-918c
Unit	30 g
Deliverydetails	No Dangerous Good /not restricted

**Description**

This Standard Reference Material (SRM) is intended for use as an analytical standard of known purity. It is intended primarily for use in the calibration and standardization of procedures for potassium (K) and chloride (Cl<sup>-</sup>) determinations employed in clinical analysis, including those using ion-selective electrodes, and for routine critical evaluation of the daily working standards used in these procedures.