

### **Tool steel C85 W2, no. 1.1630 - chips**

Art. ID                      ECRM-D 036-1  
Unit                         100 g  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	0,858	%		
	Silicon (Si)	[7440-21-3]	0,194	%		
	Manganese (Mn)	[7439-96-5]	0,327	%		
	Phosphorus (P)	[7723-14-0]	0,0074	%		
	Sulfur (S)	[7704-34-9]	0,0095	%		
	Arsenic (As)	[7440-38-2]	0,0233	%		
	Copper (Cu)	[7440-50-8]	0,065	%		
	Nitrogen (N)	[7727-37-9]	0,01	%		
	Alloy / Werkstoff		1.1630			